**PH-WERT MESSGERÄT PH100 ATC**

BEST-NR. 10 11 45

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**PH-VALUE MEASURING DEVICE PH100 ATC**

ITEM NO. 10 11 45

These operating instructions are a part of the product. They contain important notices about implementation and handling. Please bear this in mind in case you pass on the product to a third party.

Therefore, keep them for future reference!

**1. INTRODUCTION**

Dear Customer,

In purchasing this Voltcraft® product, you have made a very good decision for which we should like to thank you.

You have acquired an above-average quality product from a brand family which has distinguished itself in the field of measuring, charging and network technology by particular competence and permanent innovation.

With Voltcraft®, you will be able to tackle difficult tasks, no matter whether you are an ambitious hobbyist or a professional user. Voltcraft® offers you reliable technology at a very favourable cost-performance ratio.

We are certain: Starting to use Voltcraft will also be the commencement of a long, successful relationship.

We hope you enjoy your new Voltcraft® product!

The device was manufactured according to the latest state of technology and meets the demands of the valid European and national guidelines. The corresponding declarations and documents have been deposited with the manufacturer.

**2. PRODUCT DESCRIPTION**

The PH-100ATC is a handy pH-value measuring device with a resolution of 0.01 pH. With this measuring device, you can easily detect how acidic or alkaline a liquid is. The device has a fast response time and low power consumption. This makes it highly effective. With the hold function, you can "freeze" the measuring value on the large LCD display. The automatic temperature compensation (ATC) ensures correct measuring values even in case of temperature fluctuations. The device is suited for domestic applications as well as for (fish) ponds, swimming pools, food inspection, quality control, photo labs, schools and universities, horticulture, etc. The device is not suited for industrial application (e.g. electroplating).

**3. INTENDED USE**

Measuring the pH-value from 0.0 to 14.0 of dead, non-flammable and/or non-caustic liquids (immersion depth min. 10 resp. max. 110 mm); with automatic temperature compensation ATC.

Only a 9-V block battery may be used as voltage source.

Operation is only permitted in a dry environment; contact with moisture must be avoided at all times. Only a 9-V block battery may be used as voltage source.

Operation is not permitted under unfavourable ambient conditions. Unfavourable ambient conditions include: dust and inflammable gases, vapours or solvents, storms or stormy conditions such as strong electrostatic fields, etc.

**Use other than that described above will lead to damage to the product. In addition, this may lead to hazards such as short circuits or fire. The product as a whole may not be opened, modified or reconstructed.**

**4. SCOPE OF DELIVERY**

- PH-100ATC
- pH-electrode with storage solution
- Operating instructions

**5. SAFETY INSTRUCTIONS AND HAZARD WARNINGS**

This device left the factory in perfect condition in terms of safety engineering.

To maintain this condition and ensure safe operation, you, as the user, must comply with the safety instructions and warnings contained in these instructions.

Measuring devices do not belong in the hands of children!

In commercial institutions, the accident prevention regulations of the Employer's Liability Insurance Association for Electrical Systems and Operating Materials are to be observed.

Check the measuring probe for damage before each measuring operation!

Perception of the optical display may be impaired under adverse lighting conditions, e.g. direct sunlight.

In case of arbitrary mechanical deformation or modifications of the measuring device's electronic components, the warranty will be void.

Be careful when handling flammable or caustic liquids. wear safety gloves, goggles and an apron. Only measure in well-ventilated environments!

Only immerse the sensor head in the dead measuring media! The immersion depth may not undercut 10 mm or exceed 110 mm.

**6. INSERTING/CHANGING THE BATTERIES**

For operation, the PH-100ATC needs one battery (type 9V block, not included in the delivery).

**Proceed as follows to insert or change the batteries:**

- Unscrew the battery compartment screw on the back of the device and take off the battery compartment cover.
- Insert a new 9V block battery or replace the spent battery with a new one of the same type. If possible, use alkaline batteries. They guarantee a longer service life.
- Then, carefully close and screw down the battery compartment cover again.

It is necessary to change the batteries when the battery icon appears on the LCD. This icon indicates that the battery is spent.

Once this battery icon appears on the LC display, the battery must be changed to ensure sufficient measuring accuracy.

To prevent damage to the device from leaking batteries, remove the batteries from the device if you will not be using it for a longer period. For the same reason, we recommend that you remove flat batteries immediately.

**7. COMMISSIONING, CALIBRATION, MEASURING****Calibration**

Prior to initial operation of the PH-100ATC, we recommend calibration to adapt the pH-electrode to the measuring device. For calibration, you need a tiny screwdriver as well as calibration solutions (buffer solutions). Calibration can be performed as one-point or two-point calibration. For the more exact two-point calibration, you need two buffer solutions (pH 7.00 and pH 4.00); for one-point calibration, you only need one buffer solution. For one-point calibration, we recommend selecting a buffer solution close to the later main measuring range.

**Two-point calibration**

- Connect the pH-electrode with the BNC input socket on the top part of the PH-100ATC. Arrest the electrode plug on the measuring device with a quarter turn.
- Take the pH-electrode out of the storage solution and rinse it with distilled or de-ionised water (also after each application/measuring operation) and dry it off.
- Activate the PH-100ATC by pressing the button "ON".
- Dip the pH-electrode into a 7.00 pH-buffer/calibration solution, stir it slightly and wait until the display has stabilised (up to 5 minutes).
- By turning the calibration trimmer "CAL PH7" on the side of the PH-100ATC (below the cover) with a screwdriver, set the display value to "7.00".
- Clean the pH-electrode with distilled or de-ionised water again and wipe it dry.
- Dip the pH-electrode into a 4.00 pH-buffer/calibration solution, stir it slightly and wait until the display has stabilised (up to 5 minutes).
- By turning the calibration trimmer "SLOPE PH4/PH10" on the side of the PH-100ATC (below the cover) with a screwdriver, set the display value to "4.00".
- Repeat this step at least twice.

**One-point calibration**

- Connect the pH-electrode with the BNC input socket on the top part of the PH-100ATC. Arrest the electrode plug on the measuring device with a quarter turn.
- Take the pH-electrode out of the storage solution and rinse it with distilled or de-ionised water (also after each application/measuring operation) and dry it off.
- Activate the PH-100ATC by pressing the button "ON".
- Dip the pH-electrode into a pH-buffer/calibration solution, stir it slightly and wait until the display has stabilised (up to 5 minutes).
- By turning the calibration trimmer "CAL PH7" on the side of the PH-100ATC (below the cover) with a screwdriver, set the display value to the pH-value of the buffer solution.
- Clean the pH-electrode with distilled or de-ionised water again and wipe it dry.

**Calibration notices**

It is not necessary to calibrate the PH-100ATC prior to each measuring operation. We recommend performing a calibration prior to each tenth measuring operation or every two weeks. The buffer solution should have a similar temperature as the measuring solution.

**Measuring**

After calibration, you can perform measuring operations with the PH-100.

- Connect the pH-electrode with the BNC input socket on the top part of the PH-100ATC. Arrest the electrode plug on the measuring device with a quarter turn.
- Take the pH-electrode out of the storage solution and rinse it with distilled or de-ionised water and dry it off.
- Activate the PH-100ATC by pressing the button "ON".
- Dip the pH-electrode into the liquid to be measured and stir it briefly. Once the display has stabilized, you can read off the pH-value of the liquid directly. The automatic temperature compensation (ATC) also ensures accurate measuring values in case of different liquid temperatures. You can "freeze" the measuring value by pressing "HOLD".
- Turn the device off after measuring by pressing "OFF".
- To ensure the function of the device over a long period of time, you must clean the pH-electrode after each measuring operation (see "Maintenance and cleaning")
- Always keep the electrode in the storage solution after measuring!

#### Notice on pH-electrodes

The pH-electrode must always be kept moist to provide accurate measuring results over a long period of time. For storage, the electrode must always be kept in a storage solution. If the storage solution is used up, you can purchase a new batch (optionally). The storage solution may oxidise on the surface of the electrode when it spills. Remove these spills with distilled water. Therefore, tightly close the storage bottle when storing the electrode (do not remove the O-ring).

Air bubbles may form in the electrode during transport. In this case, shake the electrode like a fever thermometer to remove the air bubbles.

Only use distilled or de-ionised water to clean the pH-electrode. Tap water contains salt and other components, which may falsify the later measuring result or the calibration.

**The pH-electrode is a wear-and-tear part. Worn electrodes have a higher stabilisation time and a higher measuring inaccuracy. Worn electrodes are exempt from warranty. The pH-electrode is only suited for liquids.**

## 8. MAINTENANCE AND CLEANING

For cleaning (rinsing) the electrode, only use distilled or de-ionised water and dab dry with paper cloths. Do not use detergents that contain carbon, petrol, alcohol or similar substances for cleaning purposes. The surface of the measuring instrument will be corroded. The vapours are also detrimental to health and are explosive. Sharp-edged tools, screwdrivers or metal brushes may not be used for cleaning.

## 9. DISPOSAL

Old electronic devices are special waste and may not be disposed of in the household waste. When the device has become unusable, dispose of it at the community collection point in accordance with the statutory regulations in force.

## 10. TECHNICAL DATA

pH-measuring range.....	: 0.00 to 14.00 pH
Resolution.....	: 0.01 pH
Accuracy.....	: ± 0.07 pH (pH 5 – pH 9)
	: ± 0.1 pH (pH 4 – pH 10)
	: ± 0.2 pH (pH 1 – pH 4)
	: ± 0.2 pH (pH 10 – pH 14)
Measuring rate.....	: 0.4 seconds
Operating temperature .....	: 0 to max. +50°C
Working humidity .....	: < 80% r.h.
Voltage supply .....	: 9 V compound battery
Dimensions.....	: 205 x 68 x 30 mm
Weight .....	: 220g w/o electrode

#### Legal notice

These operating instructions are a publication by Voltcraft®, Lindenweg 15, D-92242 Hirschau/Germany, Phone +49 180/586 582 7 ([www.voltcraft.de](http://www.voltcraft.de)).

All rights including translation reserved. Reproduction by any method, e.g. photocopy, microfilming, or the capture in electronic data processing systems require the prior written approval by the editor. Reprinting, also in part, is prohibited.

These operating instructions represent the technical status at the time of printing. Changes in technology and equipment reserved.

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© Copyright 2012 by Voltcraft®
V3\_1112\_01/AB